



ALL TERRAIN WHEELCHAIR MODEL:AY-YSB1.0



**Function
3-in-1**



**Waterproof
material**



Recycling



**Inflatable
wheels**




**Breathable
seat back**



**120KG
max weight**



**Beach
Free travel**



DEFINITION DU PRODUIT

The All-Terrain 4-Wheel Beach Wheelchair is built for easy mobility on soft and uneven surfaces like sand, grass, and gravel. With its four wide wheels and rust-resistant frame, it provides a smooth, stable ride outdoors, helping users enjoy nature with comfort and ease.

TECHNICAL CHARACTERISTICS

Material:	Carbon steel + epoxy paint frame, stainless steel parts
Overall Size:	120 × 110 × 122 cm
Product Length:	120 cm
Seat width:	45 cm
Features:	Removable push-rest, high backrest
Front wheel:	9" PU
Rear wheel:	12" PU
Seat Fabric	Hypoallergenic, UV-resistant, waterproof, ventilated
Frame:	Folding carbon steel frame, stainless steel fork
Weight Capacity	120 kg
Backrest	90°-180° backward tilt, maintain any angle
Leg support	5-level adjustable leg rest
Seat Color	Blue
Package Size	81 × 27 × 93 cm
Net Weight:	19.5 kg
Gross Weight	24 kg

PRODUCT PARAMETERS



Be safe in the water

Maximum load-bearing capacity: 100KG
Product weight: 40KG

Length: 120cm

Height: 122cm

Self-locking brake

Removable high backrest

Seat width: 40 cm
Seat depth: 43 cm
Seat height: 49 cm

Folding backrest

5-level adjustable leg rest

Backrest 90° -180° for lying down







5-level adjustable leg rest

PU pneumatic tire: 12 inches








PU pneumatic tire: 9 inches








Note: Do not enter the water alone

PRODUCT DETAILS



INSTALLATION GUIDE







FOSHAN ARY MEDICAL CO.,LTD



Kitty Tong

Supply Chain Director

 +8613928274315

 arymedical1688@gmail.com

 Room 412, Junfeng Building, No. 9,
Tiyu 3rd Street, Zumiao Street, Foshan
City, Guangdong PR, China. Zip 528000

 <https://arymd.com>
<https://arymedical.com>
<https://ary1688.en.alibaba.com>