



WATER WHEELCHAIR MODEL:AY-ST-V1.0



**Function
3-in-1**



**Waterproof
material**



Recycling



**Inflatable
wheels**




**Breathable
seat back**



**150KG
max weight**



**Beach
Free travel**



DEFINITION DU PRODUIT

The Waterwheels is an amphibious wheelchair that makes it easier for people with reduced mobility to move around on soft ground such as sand. Thanks to its floating armrests and three wheels, it allows you to float safely in calm waters and enjoy the joys of water.

TECHNICAL CHARACTERISTICS


Material:	Aluminum alloy frame
Handrail:	PVC Inflatable Handrail
Seat cushion back:	PVC breathable UV protection
Product size:	210*110*117cm
Folding size:	145*40*76cm
Seat width:	46cm
Seat depth:	46cm
Seat height:	40-46cm
Armrest width:	60cm
Product weight:	27kg
Front wheel:	12"*4.5" - PVC pneumatic
Rear wheel:	16"*6.5" - PVC pneumatic
Backrest:	Three-level adjustment, 110°-130°
Weight capacity:	120kg (300lbs)
Packing method:	1pcs/carton
Packing size:	carton 150*40*76cm
GW KG:	36.5kg
color:	Yellow, blue, white

PRODUCT PARAMETERS











Be safe in the water

Maximum load-bearing capacity: 150KG
Product weight: 30KG

















NOTE: UV and water resistant fabric - seawater resistant material

PRODUCT DETAILS



INSTALLATION GUIDE







FOSHAN ARY MEDICAL CO.,LTD



Kitty Tong
Supply Chain Director

 +8613928274315

 arymedical1688@gmail.com

 Room 412, Junfeng Building, No. 9,
Tiyu 3rd Street, Zumiao Street, Foshan
City, Guangdong PR, China. Zip 528000

 <https://arymd.com>
<https://arymedical.com>
<https://ary1688.en.alibaba.com>