



ALL TERRAIN WHEELCHAIR MODEL:AY-ST-V3.0



**Function
3-in-1**



**Waterproof
material**



Recycling



**Inflatable
wheels**




**Breathable
seat back**



**100KG
max weight**



**Beach
Free travel**



DEFINITION DU PRODUIT

The Lightweight Beach Wheelchair is designed for effortless mobility on sandy beaches, grassy parks, and uneven outdoor terrain. Featuring wide balloon wheels and a sturdy yet lightweight frame, it ensures a smooth and stable ride while being easy to transport and maneuver. This wheelchair allows users to enjoy seaside activities and outdoor leisure with comfort, safety, and independence.

TECHNICAL CHARACTERISTICS

Material:	Carbon steel + epoxy paint frame, stainless steel parts
Overall Size:	78 x 80 x 95 cm
Product Length:	92 cm
Width Between Armrests:	42 cm
Width Between Rear Wheels:	50 cm
Product Height:	95 cm
Features:	Foldable, lightweight
Front wheel:	9"
Rear wheel:	12"
Seat Fabric	UV-resistant, waterproof, breathable
Frame:	Foldable carbon steel frame, stainless steel front fork
Weight Capacity	100 kg
Leg support	Flip-up pedals
Seat Color	Blue-black
Package Size	58 x 36 x 83 cm
Net Weight:	12.6 kg
Gross Weight	15.1 kg

PRODUCT PARAMETERS



Be safe in the water

Maximum load-bearing capacity: 100KG
Product weight: 13KG








Length: 92cm

Height: 95cm
















Note: Do not enter the water alone

PRODUCT DETAILS



INSTALLATION GUIDE





MAINTAINANCE





FOSHAN ARY MEDICAL CO.,LTD



Kitty Tong
Supply Chain Director

 +8613928274315

 arymedical1688@gmail.com

 Room 412, Junfeng Building, No. 9,
Tiyu 3rd Street, Zumiao Street, Foshan
City, Guangdong PR, China. Zip 528000

 <https://arymd.com>
<https://arymedical.com>
<https://ary1688.en.alibaba.com>