



BEACH WALKING AID MODEL:AY-STZ1.0



**Function
3-in-1**



**Waterproof
material**



Recycling



**Inflatable
wheels**



Comfort seat



**100KG
max weight**



**Beach
Free travel**



DEFINITION DU PRODUIT

The Beach Walker is designed for easy and stable movement on sandy beaches, grassy parks, and uneven outdoor terrain. Featuring wide balloon wheels and a durable lightweight frame, it ensures smooth navigation without sinking into the sand. This walker allows users to enjoy beachside strolls and outdoor leisure with comfort, safety, and confidence, offering reliable support for every step of their coastal adventure.

TECHNICAL CHARACTERISTICS

Material:	Carbon steel with epoxy paint frame
Overall Size:	95-105 x 80 x 80 cm
Product Length:	80 cm
Seat height:	55cm
Width Between Armrests:	45 cm
Width Between Rear Wheels:	50 cm
Product Height:	95-105 cm
Features:	Foldable, lightweight, quick-release version
Armrest	Height adjustable, with brake
Front wheel:	9"
Rear wheel:	12"
Seat Fabric	PU comfort seat
Frame:	Foldable carbon steel frame, stainless steel front fork
Weight Capacity	100 kg
Seat Color	black
Package Size	53 x 28 x 79 cm
Net Weight:	13.5 kg
Gross Weight	15.6 kg

PRODUCT PARAMETERS



Be safe in the water

Maximum load-bearing capacity: 100KG
Product weight: 13.5KG








Length: 80cm

Height: 95-105cm



Note: Do not enter the water alone

PRODUCT DETAILS



INSTALLATION GUIDE

